



Animal and Fungi DNA Preparation - Solution Kit

Solution based genomic DNA purification from animal tissue and fungi

Cat. No.	Amount
PP-208S	100 preparations
PP-208L	5x 100 preparations

For general laboratory use.

Shipping: shipped at ambient temperature

Storage Conditions: store at ambient temperature

Shelf Life: 12 months

Description:

Animal and Fungi DNA Preparation Kit is designed for convenient and fast isolation of genomic DNA from animal tissue and fungi. The solution based system minimizes DNA fragmentation that may be problematic in spin-column / filtration based methods. Because phenol or chloroform is not used it is safe and does not produce any harmful waste.

Solution based genomic DNA purification kits guarantee minimal DNA fragmentation and yield DNA sized up to 150 kb.

Content:

Cell Lysis Solution

Cell Resuspension Solution

Proteinase K (before use, solve in double distilled water to obtain a final concentration of 20 mg/ml) - store at -20 °C

Protein Precipitation Solution

Washing Buffer (before use, add 96-99 % Ethanol as indicated on the bottle)

DNA Hydration Solution

RNase A (before use, solve in double distilled water to obtain a final concentration of 4 mg/ml) - store at -20 °C

To be provided by you:

Isopropanol (2-propanol) >99 %

96-99 % Ethanol

Microtubes 1.5 ml

Preparation procedure:

Before start, provide >99 % Isopropanol (2-propanol) (not included in the kit).

For S pack (100 preps): Add 200 µl dd-water to the Proteinase K tube, 200 µl dd-water to the RNase A tube and 48 ml 96-99 % Ethanol (not included in the kit) to the Washing Buffer bottle.



Animal and Fungi DNA Preparation - Solution Kit

Solution based genomic DNA purification from animal tissue and fungi

Buffer	PP-208S 100 preps	PP-208L 5x 100 preps
Cell Lysis Solution	32 ml	5x 32 ml
Cell Resuspension Solution	32 ml	5x 32 ml
Proteinase K (20 mg/ml)	4 mg	5x 4 mg
Protein Precipitation Solution	11 ml	5x 11 ml
Washing Buffer	add 48 ml Ethanol (final volume 60 ml)	add 48 ml Ethanol to each bottle (final volume 60 ml each)
DNA Hydration Solution	11 ml	5x 11 ml
RNase A (4 mg/ml)	0.8 mg	5x 0.8 mg

pellet that ranges in color from off-white to light green.

- Discard the supernatant and drain tube briefly on clean absorbent paper.
- Add 500 µl Washing Buffer and invert tube several times to wash the DNA pellet.
- Centrifuge at 15,000 g for 1 min. Discard the ethanol carefully.
- Air dry at room temperature for 10-15 min.

4 DNA Hydration:

- Add 50-100 µl of DNA Hydration Solution to the dried DNA pellet.
- Add 1.5 µl of RNase A to the cell lysate.
- Mix the sample by inverting the tube and incubate at 37 °C for 30-60 min.
- Hydrate the DNA by incubating sample at 65 °C for 60 min.
- Store DNA at 4 °C. For long time storage, place sample at -20 °C or -80 °C.

1a Cell Lysis for Animal Tissue:

- Transfer 5-10 mg of fresh or frozen tissue to a 1.5 ml microtube.
- Add 300 µl Cell Lysis Solution to the tissue.
- Add 1.5 µl Proteinase K Solution to the lysate and mix by inverting several times.
- Incubate at 55 °C overnight or until tissue has dissolved.

1b Cell Lysis for Fungi:

- Transfer 1 ml of the cultured cells to a 1.5 ml microtube.
- Harvest the cells by centrifuging at 15,000 g for 1 min and discard supernatant.
- Resuspend the cell pellet in 300 µl Cell Resuspension Solution.
- Add 1.5 µl Proteinase K Solution and mix by inverting several times.
- Incubate at 55 °C for 60 min.
- Centrifuging at 15,000 g for 1 min and discard supernatant.
- Resuspend the pellet in 300 µl Cell Lysis Solution.

2 Protein Precipitation:

- Add 100 µl of Protein Precipitation Solution to the cell lysate.
- Mix the solution well by vortexing for 20 sec.
- Centrifuge at 15,000 g for 3 min. (The precipitated protein will be a tight pellet. If the pellet is not tight, repeat mixing, incubate on ice for 10 minutes, and then centrifuge again.)

3 DNA Precipitation:

- Transfer the supernatant to a clean 1.5 ml microtube containing 300 µl of Isopropanol >99 %.
- Mix the sample by inverting gently 50 times.
- Centrifuge at 15,000 g for 1 min. The DNA will be visible as a