



Miltenyi Biotec

REAscreen™ MAX Kit, human, FFPE, version 02 Panel composition – 146 antibody conjugates

130-129-734

For more information about single antibodies click on the respective clone.

Antigen	Clone	Fluorochrome	Antibody type	Plate	Plate position
Actin (Smooth Muscle)	REAL650	FITC	REAdye_lease™	Plate 1/2	D4
AKT Pan (PKB)	REA676	PE	REAffinity™	Plate 2/2	B7
Amyloid A	REA931	FITC	REAffinity	Plate 1/2	B8
AN2	REA989	PE	REAffinity	Plate 1/2	G11
Annexin I	REA1122	FITC	REAffinity	Plate 2/2	A7
BATF	REA486	FITC	REAffinity	Plate 1/2	D8
Bcl-10	REA1140	PE	REAffinity	Plate 2/2	A3
Bcl-2	REA872	FITC	REAffinity	Plate 1/2	B10
BMI-1	REA438	FITC	REAffinity	Plate 2/2	C6
CA-125	REAL909	PE	REAdye_lease	Plate 1/2	A7
Caldesmon	REA1120	PE	REAffinity	Plate 2/2	E1
Calponin	REAL635	PE	REAdye_lease	Plate 1/2	B4
Catalase	REAL553	PE	REAdye_lease	Plate 1/2	E4
CD100	REA316	FITC	REAffinity	Plate 2/2	C7
CD104 (Integrin β4)	REA236	FITC	REAffinity	Plate 2/2	A5
CD107a (LAMP-1)	REAL653	APC	REAdye_lease	Plate 1/2	G6
CD107b	REA1073	FITC	REAffinity	Plate 2/2	A1
CD117	REA787	PE	REAffinity	Plate 1/2	D9
CD123	REA918	FITC	REAffinity	Plate 1/2	E9
CD13	REAL771	PE	REAdye_lease	Plate 1/2	D6
CD134 (OX40)	ACT35	PE	Hybridoma	Plate 1/2	D7
CD138	REA929	FITC	REAffinity	Plate 1/2	D11
CD147	REA282	FITC	REAffinity	Plate 2/2	A6
CD15	VIMC6	FITC	Hybridoma	Plate 1/2	E7
CD158i (KIR2DS4)	JJC11.6	PE	Hybridoma	Plate 1/2	G7
CD162	REA319	FITC	REAffinity	Plate 2/2	H7
CD181 (CXCR1)	REA958	FITC	REAffinity	Plate 1/2	H11
CD182 (CXCR2)	REA208	PE	REAffinity	Plate 2/2	G7

Antigen	Clone	Fluorochrome	Antibody type	Plate	Plate position
CD183 (CXCR3)	REAL756	FITC	REAdye_lease	Plate 1/2	E6
CD195 (CCR5)	REA245	PE	REAffinity	Plate 1/2	A12
CD196 (CCR6)	REA190	PE	REAffinity	Plate 2/2	E4
CD2	REA1130	FITC	REAffinity	Plate 2/2	B2
CD20 Cytoplasmic	REAL649	APC	REAdye_lease	Plate 1/2	C4
CD202b (TIE-2)	REA198	APC	REAffinity	Plate 2/2	E7
CD204	REA460	PE	REAffinity	Plate 2/2	B6
CD209 (DC-SIGN)	REAL690	APC	REAdye_lease	Plate 1/2	C5
CD223	REA351	PE	REAffinity	Plate 2/2	F4
CD226 (DNAM-1)	REA1040	FITC	REAffinity	Plate 1/2	B12
CD233	REA368	FITC	REAffinity	Plate 2/2	D3
CD234 (DARC)	REA376	PE	REAffinity	Plate 2/2	H6
CD235a (Glycophorin A)	REAL490	FITC	REAdye_lease	Plate 1/2	F3
CD239 (BCAM)	REA276	PE	REAffinity	Plate 1/2	E11
CD240DCE	REA327	PE	REAffinity	Plate 1/2	E8
CD243 (ABCB1)	REA495	PE	REAffinity	Plate 2/2	E6
CD268	REA1115	FITC	REAffinity	Plate 2/2	C1
CD271 (LNGFR)	REAL709	APC	REAl ease®	Plate 1/2	A5
CD279 (PD1)	REA1165	PE	REAffinity	Plate 2/2	C4
CD295 (LEPR)	REA361	APC	REAffinity	Plate 2/2	F6
CD2AP	REA1147	PE	REAffinity	Plate 2/2	C3
CD3	REA1151	APC	REAffinity	Plate 1/2	B2
CD305 (LAIR-1)	REA447	PE	REAffinity	Plate 2/2	G1
CD317 (BST2)	REA202	FITC	REAffinity	Plate 2/2	F7
CD334 (FGFR4)	4FR6D3	FITC	Hybridoma	Plate 1/2	C7
CD43	REA833	FITC	REAffinity	Plate 1/2	G9
CD44	REA690	PE	REAffinity	Plate 1/2	B9
CD45	5B1	PE	Hybridoma	Plate 1/2	A1
CD45RA	REAL164	APC	REAl ease	Plate 1/2	G1
CD45RB	REA119	FITC	REAffinity	Plate 1/2	C3
CD45RO	REA611	PE	REAffinity	Plate 1/2	D1
CD5	REAL760	PE	REAdye_lease	Plate 1/2	B5
CD52	REA164	FITC	REAffinity	Plate 2/2	E5
CD57	REA769	PE	REAffinity	Plate 1/2	F9
CD65	REA1128	PE	REAffinity	Plate 2/2	A2
CD66 (CEA)	REA1158	FITC	REAffinity	Plate 2/2	D4
CD66abce	REA1230	PE	REAffinity	Plate 2/2	D6
CD66acde	REA428	PE	REAffinity	Plate 2/2	G6
CD66ace	REA889	FITC	REAffinity	Plate 1/2	E10
CD66b	REA306	PE	REAffinity	Plate 2/2	B5
CD66c	REA414	FITC	REAffinity	Plate 2/2	B4
CD66e	REA876	PE	REAffinity	Plate 1/2	C10
CD79a	REA1168	PE	REAffinity	Plate 2/2	H4
CD88 (C5AR)	REA1213	PE	REAffinity	Plate 2/2	D5
CD8a	REA1024	FITC	REAffinity	Plate 1/2	A2
CD90	REAL677	PE	REAdye_lease	Plate 1/2	F2
CD99	REAL687	FITC	REAl ease	Plate 1/2	H4
CLA	REA1101	PE	REAffinity	Plate 2/2	B1
Collagen I	REAL958	APC	REAdye_lease	Plate 1/2	B7
Cytokeratin	REAL648	APC	REAdye_lease	Plate 1/2	D2
Cytokeratin 10	REA933	FITC	REAffinity	Plate 1/2	A11
Cytokeratin 10/13	REA1138	FITC	REAffinity	Plate 2/2	H2

Antigen	Clone	Fluorochrome	Antibody type	Plate	Plate position
Cytokeratin 14	REAL831	FITC	REAdye_lease	Plate 1/2	F5
Cytokeratin 19	REAL822	PE	REAdye_lease	Plate 1/2	E5
Cytokeratin 20	REAL753	APC	REAdye_lease	Plate 1/2	H5
Cytokeratin 5/6	REAL538	PE	REAdye_lease	Plate 1/2	G3
Cytokeratin 5/6/8/17/19	REA1141	FITC	REAffinity	Plate 2/2	E3
Cytokeratin 7	REA935	FITC	REAffinity	Plate 1/2	B11
Cytokeratin 8	REA885	FITC	REAffinity	Plate 1/2	D10
Cytokeratin HMW	REAL645	FITC	REAdye_lease	Plate 1/2	C6
Dectin-1	REA515	PE	REAffinity	Plate 2/2	H5
Desmin	REA1134	FITC	REAffinity	Plate 2/2	B8
Dnmt3b	REA565	APC	REAffinity	Plate 2/2	A8
EZH2	REA907	PE	REAffinity	Plate 1/2	F10
FcεR1α	REA758	FITC	REAffinity	Plate 1/2	C9
Galectin-3	REA1072	PE	REAffinity	Plate 1/2	G12
Galectin-9	REA435	FITC	REAffinity	Plate 2/2	C5
GCDFP-15	REAL828	FITC	REAdye_lease	Plate 1/2	F4
GFAP	REAL735	PE	REAdye_lease	Plate 1/2	G5
Hepatocyte	REA1088	FITC	REAffinity	Plate 1/2	A9
HLA-ABC	REAL597	APC	REAdye_lease	Plate 1/2	A3
HLA-DR	REAL550	FITC	REAdye_lease	Plate 1/2	C2
HLA-DR, DP, DQ	REA332	FITC	REAffinity	Plate 1/2	G2
Hsp70	REA349	PE	REAffinity	Plate 1/2	C12
IgD	REA740	PE	REAffinity	Plate 1/2	B1
IRF-7	REA521	FITC	REAffinity	Plate 1/2	G8
IRF-7 pS477/pS479	REA310	PE	REAffinity	Plate 1/2	C8
Jak1	REA700	PE	REAffinity	Plate 1/2	E12
JNK2 (MAPK9)	REA1153	FITC	REAffinity	Plate 2/2	H3
Ki-67	REA1123	FITC	REAffinity	Plate 2/2	H1
LRP-4	REA552	APC	REAffinity	Plate 2/2	E2
MART-1	REA1093	FITC	REAffinity	Plate 1/2	H12
Mast Cell Tryptase	REAL798	FITC	REAdye_lease	Plate 1/2	D5
Melanocyte PMEL	REA927	PE	REAffinity	Plate 1/2	C11
Metallothionein	REA1139	FITC	REAffinity	Plate 2/2	B3
MLC2v	REA401	PE	REAffinity	Plate 2/2	G2
Myogenin	REA930	FITC	REAffinity	Plate 1/2	F11
Myosin Smooth Muscle	REA1107	FITC	REAffinity	Plate 2/2	D1
NeuN	REA1131	PE	REAffinity	Plate 2/2	C2
p120 Catenin pS879	REA548	PE	REAffinity	Plate 1/2	F1
p53	REA1132	FITC	REAffinity	Plate 2/2	D2
PAX-5 (BSAP)	REA140	PE	REAffinity	Plate 1/2	A8
PCNA	REA858	FITC	REAffinity	Plate 1/2	A10
PI3K p85 α	REAL589	PE	REAdye_lease	Plate 1/2	G4
PKCa	REA156	APC	REAffinity	Plate 2/2	C8
PLAP	REA1089	FITC	REAffinity	Plate 1/2	D12
Plasma Cell	REA908	FITC	REAffinity	Plate 1/2	G10
PLP	REA1155	PE	REAffinity	Plate 2/2	A4
Podoplanin	REAL468	PE	REAllease	Plate 1/2	E3
Slug	REA404	PE	REAffinity	Plate 1/2	H8
Sox2	REA320	FITC	REAffinity	Plate 2/2	G4
SSRP1	REA555	PE	REAffinity	Plate 1/2	F8
SUSD2	REA795	PE	REAffinity	Plate 1/2	H10
Syk	REA111	PE	REAffinity	Plate 2/2	D7

Antigen	Clone	Fluorochrome	Antibody type	Plate	Plate position
Synaptophysin	REA1121	FITC	REAffinity	Plate 2/2	F1
TIM-3	REAL818	APC	REAdye_lease	Plate 1/2	F7
TOM22	REAL439	APC	REAdye_lease	Plate 1/2	A4
TTF-1	REA1090	FITC	REAffinity	Plate 1/2	F12
Tyrosinase	REA1149	PE	REAffinity	Plate 2/2	F3
Vimentin	REAL1008	APC	REAdye_lease	Plate 1/2	D3
WT1	REAL748	APC	REAdye_lease	Plate 1/2	B6
ZAP70	REA814	PE	REAffinity	Plate 1/2	H9
α Tubulin	REA1136	FITC	REAffinity	Plate 2/2	F2
α -Actinin (Sarcomeric)	REA402	PE	REAffinity	Plate 2/2	F5
β -Actin	REAL571	PE	REAdye_lease	Plate 1/2	H3
β -Catenin	REA480	FITC	REAffinity	Plate 2/2	G5
β -Tubulin 3	REA1152	PE	REAffinity	Plate 2/2	G3
γ -Actin 1	REAL591	FITC	REAdye_lease	Plate 1/2	H6



Miltenyi Biotec B.V. & Co. KG | Phone +49 2204 8306-0 | Fax +49 2204 85197 | macsde@miltenyi.com | www.miltenyibiotec.com

Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS, the Miltenyi Biotec logo, REAdye_lease, REAffinity, REAlease, and REAscreen are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. Copyright © 2022 Miltenyi Biotec and/or its affiliates. All rights reserved.